

Time to fly

Hear what you want to hear by amplifying the sounds and voices most important to you.

- 1) Plug your earphone or headphone into the headphone jack.
- 2) Position the internal microphone near the preferred sound.
- 3) Adjust the volume to your comfort.



Take Flight

Transform your senses. Listen to life, to the sounds you choose, with the new Pocketalker® 2.0 personal amplifier from Williams Sound.

Designed with you in mind, the 2.0 exudes beauty and grace both inside and out. Clean lines designed for comfort and style. Cutting-edge technology clearly in tune with your life.

So, spread those wings. Personally tailor your listening experiences with the 2.0's superior sound and unparalleled frequency range.

Amplify your life with Pocketalker 2.0

Wherever you land



One-on-one discussions



Television listening



Conversing in the car



Restaurants

PERSONAL AMPLIFIER

1.800.843.3544 (USA)
+1.952.943.2252 (World)

info@williamssound.com (USA)

info-intl@williamssound.com (World)

www.williamssound.com

 WILLIAMS SOUND®

©Copyright 2015 - Williams Sound, LLC BROM 052A





Beauty is **more** than skin deep.



Ergonomic, lightweight design



Extended high and low frequency range



Superior sound quality



Internal microphone



Built-in T-coil loop compatible



Rechargeable ready



Hearing aids vs. **Pocketalker**

A hearing aid is a customized medical device that amplifies **all** sounds.

The Pocketalker personal amplifier is used to amplify **specific** sounds or voices. It helps you hear what you want to hear by reducing distracting background noise.



HELPING PEOPLE HEAR