


POCKETALKER COMPARISON



	PRO	Ultra	2.0
	● ● ● ● ●	● ● ● ● ●	● ● ● ● ●
Product Warranty (5-Year)	•	•	•
AAA (2) Alkaline Batteries		•	•
AA (2) Alkaline Batteries	•		
Battery Life (105 Hrs Alkaline)			•
Battery Life (200 Hrs Alkaline)	•	•	
Battery Life (40 Hrs Rechargeable)	•		
Battery Life (80 Hrs Rechargeable)			•
AAA (2) Rechargeable Batteries			•
AA (2) Rechargeable Batteries	•		
Rechargeable Ready	•		•
External Microphone	•	•	•
Internal Microphone			•
Volume Control	•	•	•
Tone Control (External)		•	•
Tone Control (Internal)	•		
Balance Adjustment			•
Built-In T-Coil			•
Mono Ear/Headphones/Neckloop	•	•	
Stereo Ear/Headphones/Neckloop			•
Includes Lanyard / Belt Clip Case	•	•	•
Shielded from Handling Noise	•	•	•
CE Approval		•	<i>Pending</i>
FDA Registered	•	•	•